



Tikrit University  
College of Veterinary Medicine

## Lecture Title: Determination of Ash in Feedstuffs

Subject name: Animal Nutrition

Subject year: Second Year

Lecturer name: Dr. Thamer Ahmed –

Dr. Ali Qaeas

Academic Email: [thamer.a.k@tu.edu.iq](mailto:thamer.a.k@tu.edu.iq)

Please do not change font type and

## Determination of Ash in Feedstuffs

### Definition:

Ash is defined as the inorganic residue that remains after burning a feedstuff sample at **600 °C** in a muffle furnace.

During incineration in the muffle furnace, all **organic matter** of the feed sample is completely oxidized, leaving behind inorganic minerals (ash). However, ash may sometimes contain **non-feed contaminants** such as stones, gravel, glass pieces, dust, debris, or sand (silica).

**Note:** In concentrate feed samples, the percentage of ash should not exceed **10%**; otherwise, the sample may be considered adulterated.

### Apparatus

- Porcelain or silica crucibles (preferred over metal cans due to their high heat tolerance).
- Muffle furnace (capable of maintaining 500–600 °C).
- Desiccator.
- Analytical balance (digital scale).

### Experimental Procedure

#### 1. Preparation of the crucible

- Place a clean, labeled crucible in the muffle furnace at **500–600 °C** for one hour.
- Transfer it to a desiccator for cooling.
- This step sterilizes the crucible and removes any adhering debris.

#### 2. Weighing before ignition

- Weigh the cooled crucible on a digital balance and record the weight (**W1 = weight of empty crucible**).
- Add **0.5 g** of finely ground feed sample. Record the combined weight (**W2 = weight before burning**).

#### 3. Ashing process

- Place the crucible with the sample in the muffle furnace at **500–600 °C** for three hours.

- The temperature must be raised gradually to avoid **spattering**, which is the ejection of sample particles caused by sudden heating.

#### 4. Cooling after ignition

- After incineration, partially cool the crucible outside the furnace, then place it in a desiccator until completely cooled.
- Weigh the crucible with ash and record (**W3 = weight after burning**).
- The crucible containing the ash can be preserved for silica determination.

**Note:** If the ash percentage exceeds **15%**, the feed sample is regarded as adulterated.

### Calculations

1. Weight of sample ( $W_s$ ) =  $W_2 - W_1$
2. Weight of ash ( $W_a$ ) =  $W_3 - W_1$
3. Percentage of ash =  $\frac{W_a}{W_s} \times 100$
4. Percentage of organic matter = Percentage of dry matter – Percentage of ash

### Example

A **2 g** sample of ground corn was incinerated in the muffle furnace at **500 °C** for three hours. After burning, the weight of ash was **0.1 g**. The percentage of dry matter of the sample is **95%**.

### Solution:

$$\text{Percentage of ash} = \frac{W_a}{W_s} \times 100$$

$$\frac{0.1}{2} \times 100 = 5\%$$

$$\text{Percentage of organic matter} = 95 - 5 = 90 \%$$

### Flowchart: Determination of Ash in Feedstuffs

#### Step 1: Preparation of Crucible

Clean crucible → Heat in muffle furnace at 500–600 °C for 1 hr → Cool in desiccator → Weigh ( $W_1$ ).

### **Step 2: Sample Weighing**

Add 0.5 g of feed sample → Weigh crucible + sample (W2).

### **Step 3: Ignition (Ashing Process)**

Place the crucible in the muffle furnace (500–600 °C) for 3 hrs. Gradual heating to prevent spattering.

### **Step 4: Cooling & Final Weighing**

Partially cool crucible → Transfer to desiccator → Weigh crucible + ash (W3).

### **Step 5: Calculations**

- Weight of sample (Ws) = W2 – W1
- Weight of ash (Wa) = W3 – W1
- % Ash = (Wa / Ws) × 100
- % Organic matter = % Dry matter – % Ash

### **Note:**

- In concentrates, ash should not exceed **10%**.
- If ash > **15%**, sample may be considered adulterated.